




Accesori Optionale pentru Distilatoare GFL

	Dechlorite Filter - GFL-2904	
	<ul style="list-style-type: none"> ○ Filtru de dechlorinare (pentru toate modele) <p>Dechlorite Filter (for all GFL Water Stills). Eliminates chlorine particles in tap water, water added by the local waterworks. Complete with connection for pressure hose ½" (inch), with first filling. (Hoses for water in- and outlet are not included)</p>	<p>Pret: - Euro F.TVA</p>

	Phosphate Cartridge - GFL-2906	
	<ul style="list-style-type: none"> ○ Filtru de fosfatare (pentru toate modele) <p>Phosphate Cartridge (for all GFL Water Stills). Prevents scale formation in the condenser by phsophatising of tap water. Complete with connection for pressure hose ½" (inch), including first filling. (Hoses for water in- and outlet are not included).</p>	<p>Pret: - Euro F.TVA</p>

	Spare filling for dechlorite filter - GFL-2905	
	<ul style="list-style-type: none"> ○ Material de umplere pentru filtrul de dechlorinare <p>Spare filling for dechlorite filter</p>	<p>Pret: - Euro F.TVA</p>

	Spare filling for phosphate cartridge - GFL-2907	
	<ul style="list-style-type: none"> ○ Material de umplere pentru filtrul de fosfatare <p>Spare filling for phosphate cartridge</p>	<p>Pret: - Euro F.TVA</p>

Separate Water Supply - GFL-2901	
<ul style="list-style-type: none"> ○ Elemente de filtrare suplimentare pentru distilatoarele din gama 2002 – 2012 pentru cresterea eficientei procesului de distilare. <p>Separate Water Supply (only for Stills 2002-2012) to feed evaporator with softened or desalinated water, and cooling coil with phosphate treated or tap water. Efficiency of the still is reduced by aprox.10-15%.</p>	<p>Pret: - Euro F.TVA</p>

Distribuitor autorizat GFL - Romania

Multi-Lab SRL
Tel. 07240 44800
Fax. +4 (021) 323 3331