


<b>Codice</b> Code Code Kode	<b>Modello</b> Model Modèle Modell	<b>Temp.°C</b> Temp.°C Temp.°C Temp.°C	<b>Capacità</b> Capacity Capacité Kapazität	<b>Prezzo €</b> Price € Prix € Preis €
---------------------------------------	---	---	--	---




AHB204

**ULF 120 W PRO 2**

-60° -86°C

120

**12.138**




AHC404

**ULF 240 W PRO 2**

-60° -86°C

240

**13.410**




AHE404

**ULF 440 W PRO 2**

-60° -86°C

440

**14.840**




AHF504

**ULF 550 W PRO 2**

-60° -86°C

550

**16.130**



AHG604

**ULF 660 W PRO 2**

-60° -86°C

660

**17.399**





- STANDARD EQUIPMENT - DOTATION - AUSSTATTUNG

Mod.	120	240	440	550	660
	75 x 79 x 133	75 x 79 x 189	90 x 95 x 193	105 x 95 x 193	118 x 95 x 193
	W	W	W	W	W
	S	S	S	S	S
	130	130	130	130	130
	4	4	4	4	4
	1	1	1	1	1
	1	3	3	3	3
	•	•	•	•	•
	•	•	•	•	•
	•	•	•	•	•
<b>GAS</b>	R404a - R508b R170	R404a - R508b R170	R404a - R508b R170	R404a - R508b R170	R404a - R508b R170
	•	•	•	•	•
	230/1/50	230/1/50	230/1/50	230/1/50	230/1/50

## ACCESSORI - ACCESSORIES - ACCESSOIRES - ZUBEHÖRTEILE

	Mod.	120	240	440	550	660
	ZRA10	€ 225.00	€ 225.00	-	-	-
	ZRA20	-	-	€ 255.00	€ 255.00	€ 255.00
	ZSC10	€ 12.00	€ 12.00	€ 12.00	€ 12.00	€ 12.00
 <b>CO2</b>	ZCO10	€ 1.930	€ 1.930	€ 1.930	€ 1.930	€ 1.930
	ZRG50	€ 2.240	€ 2.240	€ 2.240	€ 2.240	€ 2.240
	ZGU10	€ 390.00	€ 390.00	€ 390.00	€ 390.00	€ 390.00
 <b>240V / 60 Hz</b>	ZVO10	+4%	+4%	+4%	+4%	+4%
	ZIAG5	€ 195.00	-	-	-	-
	ZIAH5	-	€ 240.00	-	-	-
	ZIAI5	-	-	€ 260.00	-	-
	ZIAJ5	-	-	-	€ 280.00	-
	ZIAK5	-	-	-	-	€ 300.00

