

SIEVE SHAKERS



SIEVE SHAKER MLBA 200 N

Electromagnetic and digital sieve shaker
Ø Sieves : 200 / 203 mm.Ø
Maximum amplitude of vibration : 2,00 mm.
Movement: 3D

Price : 1.290,00 €
Code : MLTZBA200N



SIEVE SHAKER MLRP 200 N

Electromagnetic and digital sieve shaker of oscillant masses.
Ø Sieves : 200 / 203 mm.Ø
Maximum amplitude of vibration : 2,50 mm.

Price : 1.490,00 €
Code : MLTZRP200N



SIEVE SHAKER MLBA 300 N

Electromagnetic and digital sieve shaker of oscillant masses.
Ø Sieves : 200 / 315 mm.Ø
Maximum amplitude of vibration : 2,50 mm.


Price : 1.950,00 €
Code : MLTZBA300N



SIEVE SHAKER MLBA 400 N

Electromagnetic and digital sieve shaker for big quantities of sample.
Ø sieves: 200 / 450 mm.Ø
Movement : 3D

Price : 3.900,00 €
Code : MLTZBA400N



SIEVE SHAKER MLBVN 258

Electromecanic sieve shaker with tapping motion.

Ø sieves : 200 / 203 mm.Ø

Orbital and tapping movement.

Price : 2.650,00 €

Code : MLTZBVN258



SIEVE SHAKER MLAJ 20

Air jet sieve shaker.

Ø sieve : 200 mm.Ø

Price : 6.210,00 €

Code : MLTZAJ20



SIEVE SHAKER MLLP 075 E

Sonic sieve shaker

Ø sieves : 75 mm.Ø

Price : 5.875,00 €

Code : MLTZLP075E



VIBROSEPARATOR MLVIB 400

Electromecanic vibroseparator

Ø sieves : 400 mm.Ø

Price : 3.490,00 €

Code : MLTZVIB400

**TEST SIEVES IN STAINLESS STEEL
ACCORDING ISO 3310.1
WIRE MESH**

	60 mm.Ø x 18 mm.H		75 mm.Ø x 30 mm.H		75 mm.Ø x 60 mm.H		75 mm.Ø x 70 mm.H		100 mm.Ø x 23 mm.H		100 mm.Ø x 40 mm.H		150 mm.Ø x 50 mm.H	
Opening	M-LB code	Price	M-LB code	Price	M-LB code	Price	M-LB code	Price	M-LB code	Price	M-LB code	Price	M-LB code	Price
125,00 mm.#														
112,00 mm.#														
106,00 mm.#														
100,00 mm.#														
90,00 mm.#														
80,00 mm.#														
75,00 mm.#														
71,00 mm.#														
63,00 mm.#														
56,00 mm.#														
53,00 mm.#														
50,00 mm.#														
45,00 mm.#														
40,00 mm.#														
37,50 mm.#														
35,50mm.#														
31,50 mm.#														
28,00 mm.#														
26,50 mm.#														
25,00 mm.#														
22,40 mm.#														
20,00 mm.#														
19,00 mm.#														
18,00 mm.#														
16,00 mm.#													0515401	52.85 €
14,00 mm.#													0515201	52.85 €
13,20 mm.#													0515001	52.85 €
12,50 mm.#													0514801	52.85 €
11,20 mm.#													0514601	52.85 €
10,00 mm.#													0514401	52.85 €

TEST SIEVES IN STAINLESS STEEL
ACCORDING ISO 3310.1
WIRE MESH

Abertura	60 mm.Ø x 18 mm.H		75 mm.Ø x 30 mm.H		75 mm.Ø x 60 mm.H		75 mm.Ø x 70 mm.H		100 mm.Ø x 23 mm.H		100 mm.Ø x 40 mm.H		150 mm.Ø x 50 mm.H	
9,50 mm.#													0514201	52.85 €
9,00 mm.#													0514001	52.85 €
8,00 mm.#													0513801	52.85 €
7,10 mm.#													0513601	52.85 €
6,70 mm.#													0513401	52.85 €
6,30 mm.#													0513201	52.85 €
5,60 mm.#													0513001	52.85 €
5,00 mm.#													0512801	52.85 €
4,75 mm.#													0512601	52.85 €
4,50 mm.#													0512401	52.85 €
4,00 mm.#									0312201	46.20 €	0412201	46.20 €	0512201	52.85 €
3,55 mm.#									0312001	46.20 €	0412001	46.20 €	0512001	52.85 €
3,35 mm.#									0311801	46.20 €	0411801	46.20 €	0511801	52.85 €
3,15 mm.#	0111601	50.40 €	2711601	57.50 €	0211601	57.50 €	2911601	57.50 €	0311601	46.20 €	0411601	46.20 €	0511601	52.85 €
2,80 mm.#	0111401	50.40 €	2711401	57.50 €	0211401	57.50 €	2911401	57.50 €	0311401	46.20 €	0411401	46.20 €	0511401	52.85 €
2,50 mm.#	0111201	50.40 €	2711201	57.50 €	0211201	57.50 €	2911201	57.50 €	0311201	46.20 €	0411201	46.20 €	0511201	52.85 €
2,36 mm.#	0111001	50.40 €	2711001	57.50 €	0211001	57.50 €	2911001	57.50 €	0311001	46.20 €	0411001	46.20 €	0511001	52.85 €
2,24 mm.#	0110801	50.40 €	2710801	57.50 €	0210801	57.50 €	2910801	57.50 €	0310801	46.20 €	0410801	46.20 €	0510801	52.85 €
2,00 mm.#	0110601	50.40 €	2710601	57.50 €	0210601	57.50 €	2910601	57.50 €	0310601	46.20 €	0410601	46.20 €	0510601	52.85 €
1,80 mm.#	0110401	50.40 €	2710401	57.50 €	0210401	57.50 €	2910401	57.50 €	0310401	46.20 €	0410401	46.20 €	0510401	52.85 €
1,70 mm.#	0110201	50.40 €	2710201	57.50 €	0210201	57.50 €	2910201	57.50 €	0310201	46.20 €	0410201	46.20 €	0510201	52.85 €
1,60 mm.#	0110001	50.40 €	2710001	57.50 €	0210001	57.50 €	2910001	57.50 €	0310001	46.20 €	0410001	46.20 €	0510001	52.85 €
1,40 mm.#	0109801	50.40 €	2709801	57.50 €	0209801	57.50 €	2909801	57.50 €	0309801	46.20 €	0409801	46.20 €	0509801	52.85 €
1,25 mm.#	0109601	50.40 €	2709601	57.50 €	0209601	57.50 €	2909601	57.50 €	0309601	46.20 €	0409601	46.20 €	0509601	52.85 €
1,18 mm.#	0109401	50.40 €	2709401	57.50 €	0209401	57.50 €	2909401	57.50 €	0309401	46.20 €	0409401	46.20 €	0509401	52.85 €
1,12 mm.#	0109201	50.40 €	2709201	57.50 €	0209201	57.50 €	2909201	57.50 €	0309201	46.20 €	0409201	46.20 €	0509201	52.85 €
1,00 mm.#	0109001	50.40 €	2709001	57.50 €	0209001	57.50 €	2909001	57.50 €	0309001	46.20 €	0409001	46.20 €	0509001	52.85 €
900 µm.#	0108801	50.40 €	2708801	57.50 €	0208801	57.50 €	2908801	57.50 €	0308801	46.20 €	0408801	46.20 €	0508801	52.85 €
850 µm.#	0108601	50.40 €	2708601	57.50 €	0208601	57.50 €	2908601	57.50 €	0308601	46.20 €	0408601	46.20 €	0508601	52.85 €

**TEST SIEVES IN STAINLESS STEEL
ACCORDING ISO 3310.1
WIRE MESH**

Opening	60 mm.Ø x 18 mm.H		75 mm.Ø x 30 mm.H		75 mm.Ø x 60 mm.H		75 mm.Ø x 70 mm.H		100 mm.Ø x 23 mm.H		100 mm.Ø x 40 mm.H		150 mm.Ø x 50 mm.H	
	Code	Price	Code	Price	Code	Price	Code	Price	Code	Price	Code	Price	Code	Price
800 µm.#	0108401	50.40 €	2708401	57.50 €	0208401	57.50 €	2908401	57.50 €	0308401	46.20 €	0408401	46.20 €	0508401	52.85 €
710 µm.#	0108201	50.40 €	2708201	57.50 €	0208201	57.50 €	2908201	57.50 €	0308201	46.20 €	0408201	46.20 €	0508201	52.85 €
630 µm.#	0108001	50.40 €	2708001	57.50 €	0208001	57.50 €	2908001	57.50 €	0308001	46.20 €	0408001	46.20 €	0508001	52.85 €
600 µm.#	0107801	50.40 €	2707801	57.50 €	0207801	57.50 €	2907801	57.50 €	0307801	46.20 €	0407801	46.20 €	0507801	52.85 €
560 µm.#	0107601	50.40 €	2707601	57.50 €	0207601	57.50 €	2907601	57.50 €	0307601	46.20 €	0407601	46.20 €	0507601	52.85 €
500 µm.#	0107401	50.40 €	2707401	57.50 €	0207401	57.50 €	2907401	57.50 €	0307401	46.20 €	0407401	46.20 €	0507401	52.85 €
450 µm.#	0107201	50.40 €	2707201	57.50 €	0207201	57.50 €	2907201	57.50 €	0307201	46.20 €	0407201	46.20 €	0507201	52.85 €
425 µm.#	0107001	50.40 €	2707001	57.50 €	0207001	57.50 €	2907001	57.50 €	0307001	46.20 €	0407001	46.20 €	0507001	52.85 €
400 µm.#	0106801	50.40 €	2706801	57.50 €	0206801	57.50 €	2906801	57.50 €	0306801	46.20 €	0406801	46.20 €	0506801	52.85 €
355 µm.#	0106601	50.40 €	2706601	57.50 €	0206601	57.50 €	2906601	57.50 €	0306601	46.20 €	0406601	46.20 €	0506601	52.85 €
315 µm.#	0106401	50.40 €	2706401	57.50 €	0206401	57.50 €	2906401	57.50 €	0306401	46.20 €	0406401	46.20 €	0506401	52.85 €
300 µm.#	0106201	50.40 €	2706201	57.50 €	0206201	57.50 €	2906201	57.50 €	0306201	46.20 €	0406201	46.20 €	0506201	52.85 €
280 µm.#	0106001	50.40 €	2706001	57.50 €	0206001	57.50 €	2906001	57.50 €	0306001	46.20 €	0406001	46.20 €	0506001	52.85 €
250 µm.#	0105801	50.40 €	2705801	57.50 €	0205801	57.50 €	2905801	57.50 €	0305801	46.20 €	0405801	46.20 €	0505801	52.85 €
224 µm.#	0105601	50.40 €	2705601	57.50 €	0205601	57.50 €	2905601	57.50 €	0305601	46.20 €	0405601	46.20 €	0505601	52.85 €
212 µm.#	0105401	50.40 €	2705401	57.50 €	0205401	57.50 €	2905401	57.50 €	0305401	46.20 €	0405401	46.20 €	0505401	52.85 €
200 µm.#	0105201	50.40 €	2705201	57.50 €	0205201	57.50 €	2905201	57.50 €	0305201	46.20 €	0405201	46.20 €	0505201	52.85 €
180 µm.#	0105001	50.40 €	2705001	57.50 €	0205001	57.50 €	2905001	57.50 €	0305001	46.20 €	0405001	46.20 €	0505001	52.85 €
160 µm.#	0104801	50.40 €	2704801	57.50 €	0204801	57.50 €	2904801	57.50 €	0304801	46.20 €	0404801	46.20 €	0504801	52.85 €
150 µm.#	0104601	50.40 €	2704601	57.50 €	0204601	57.50 €	2904601	57.50 €	0304601	46.20 €	0404601	46.20 €	0504601	52.85 €
140 µm.#	0104401	50.40 €	2704401	57.50 €	0204401	57.50 €	2904401	57.50 €	0304401	46.20 €	0404401	46.20 €	0504401	52.85 €
125 µm.#	0104201	50.40 €	2704201	57.50 €	0204201	57.50 €	2904201	57.50 €	0304201	46.20 €	0404201	46.20 €	0504201	52.85 €
112 µm.#	0104001	53.00 €	2704001	57.50 €	0204001	57.50 €	2904001	57.50 €	0304001	46.20 €	0404001	46.20 €	0504001	52.85 €
106 µm.#	0103801	53.00 €	2703801	57.50 €	0203801	57.50 €	2903801	57.50 €	0303801	46.20 €	0403801	46.20 €	0503801	52.85 €
100 µm.#	0103601	53.00 €	2703601	57.50 €	0203601	57.50 €	2903601	57.50 €	0303601	46.20 €	0403601	46.20 €	0503601	52.85 €
90 µm.#	0103401	53.00 €	2703401	57.50 €	0203401	57.50 €	2903401	57.50 €	0303401	46.20 €	0403401	46.20 €	0503401	52.85 €
80 µm.#	0103201	53.00 €	2703201	57.50 €	0203201	57.50 €	2903201	57.50 €	0303201	46.20 €	0403201	46.20 €	0503201	52.85 €
75 µm.#	0103001	53.00 €	2703001	57.50 €	0203001	57.50 €	2903001	57.50 €	0303001	46.20 €	0403001	46.20 €	0503001	52.85 €
71 µm.#	0102801	53.00 €	2702801	57.50 €	0202801	57.50 €	2902801	57.50 €	0302801	46.20 €	0402801	46.20 €	0502801	52.85 €

**TEST SIEVES IN STAINLESS STEEL
ACCORDING ISO 3310.1
WIRE MESH**

Abertura	60 mm.Ø x 18 mm.H		75 mm.Ø x 30 mm.H		75 mm.Ø x 60 mm.H		75 mm.Ø x 70 mm.H		100 mm.Ø x 23 mm.H		100 mm.Ø x 40 mm.H		150 mm.Ø x 50 mm.H	
	63 µm.#	0102601	53.00 €	2702601	68.25 €	0202601	68.25 €	2902601	68.25 €	0302601	46.20 €	0402601	46.20 €	0502601
56 µm.#	0102401	53.00 €	2702401	68.25 €	0202401	68.25 €	2902401	68.25 €	0302401	52.50 €	0402401	52.50 €	0502401	64.85 €
53 µm.#	0102201	53.00 €	2702201	68.25 €	0202201	68.25 €	2902201	68.25 €	0302201	52.50 €	0402201	52.50 €	0502201	64.85 €
50 µm.#	0102001	53.00 €	2702001	68.25 €	0202001	68.25 €	2902001	68.25 €	0302001	52.50 €	0402001	52.50 €	0502001	64.85 €
45 µm.#	0101801	53.00 €	2701801	68.25 €	0201801	68.25 €	2901801	68.25 €	0301801	52.50 €	0401801	52.50 €	0501801	64.85 €
40 µm.#	0101601	53.00 €	2701601	68.25 €	0201601	68.25 €	2901601	68.25 €	0301601	62.00 €	0401601	62.00 €	0501601	71.80 €
38 µm.#	0101401	77.00 €	2701401	84.00 €	0201401	84.00 €	2901401	84.00 €	0301401	62.00 €	0401401	62.00 €	0501401	71.80 €
36 µm.#	0101201	77.00 €	2701201	84.00 €	0201201	84.00 €	2901201	84.00 €	0301201	62.00 €	0401201	62.00 €	0501201	71.80 €
32 µm.#	0101001	77.00 €	2701001	84.00 €	0201001	84.00 €	2901001	84.00 €	0301001	69.85 €	0401001	69.85 €	0501001	88.00 €
25 µm.#	0100801	92.00 €	2700801	99.75 €	0200801	99.75 €	2900801	99.75 €	0300801	74.55 €	0400801	74.55 €	0500801	125.00 €
20 µm.#	0100601	148.00 €	2700601	168.00 €	0200601	168.00 €	2900601	168.00 €	0300601	153.00 €	0400601	153.00 €	0500601	315.00 €

SPECIAL OPENINGS (NON STANDARDIZED)

Opening	60 mm.Ø x 18 mm.H		75 mm.Ø x 30 mm.H		75 mm.Ø x 60 mm.H		75 mm.Ø x 70 mm.H		100 mm.Ø x 23 mm.H		100 mm.Ø x 40 mm.H		150 mm.Ø x 50 mm.H	
	6,00 mm.#			2713101	57.50 €	0213101	57.50 €	2913101	57.50 €	0313101	46.20 €	0413101	46.20 €	0513101
3,00 mm.#	0111501	50.40 €	2711501	57.50 €	0211501	57.50 €	2911501	57.50 €	0311501	46.20 €	0411501	46.20 €	0511501	52.85 €
1,50 mm.#	0109901	50.40 €	2709901	57.50 €	0209901	57.50 €	2909901	57.50 €	0309901	46.20 €	0409901	46.20 €	0509901	52.85 €
750 µm.#	0108301	50.40 €	2708301	57.50 €	0208301	57.50 €	2908301	57.50 €	0308301	46.20 €	0408301	46.20 €	0508301	52.85 €

TEST SIEVES IN STAINLESS STEEL
ACCORDING ISO 3310.1
WIRE MESH

Opening	200 mm.Ø Air Jet 25 mm.H		200 mm.Ø New System 25 mm.H		200 mm.Ø x 25 mm.H		200 mm.Ø x 50 mm.H		203 mm.Ø x 25 mm.H		203 mm.Ø x 50 mm.H		250 mm.Ø x 70 mm.H	
	M-LB code	Price	M-LB code	Price	M-LB code	Price	M-LB code	Price	M-LB code	Price	M-LB code	Price	M-LB code	Price
125,00 mm.#							0820201	62.60 €			1020201	62.60 €	1220201	99.50 €
112,00 mm.#							0820001	62.60 €			1020001	62.60 €	1220001	99.50 €
106,00 mm.#							0819801	62.60 €			1019801	62.60 €	1219801	99.50 €
100,00 mm.#							0819601	62.60 €			1019601	62.60 €	1219601	99.50 €
90,00 mm.#							0819401	62.60 €			1019401	62.60 €	1219401	99.50 €
80,00 mm.#							0819201	62.60 €			1019201	62.60 €	1219201	99.50 €
75,00 mm.#							0819001	62.60 €			1019001	62.60 €	1219001	99.50 €
71,00 mm.#							0818801	62.60 €			1018801	62.60 €	1218801	99.50 €
63,00 mm.#							0818601	62.60 €			1018601	62.60 €	1218601	99.50 €
56,00 mm.#							0818401	62.60 €			1018401	62.60 €	1218401	99.50 €
53,00 mm.#							0818201	62.60 €			1018201	62.60 €	1218201	99.50 €
50,00 mm.#							0818001	62.60 €			1018001	62.60 €	1218001	99.50 €
45,00 mm.#							0817801	62.60 €			1017801	62.60 €	1217801	99.50 €
40,00 mm.#							0817601	62.60 €			1017601	62.60 €	1217601	99.50 €
37,50 mm.#							0817401	62.60 €			1017401	62.60 €	1217401	99.50 €
35,50mm.#							0817201	62.60 €			1017201	62.60 €	1217201	99.50 €
31,50 mm.#							0817001	62.60 €			1017001	62.60 €	1217001	99.50 €
28,00 mm.#							0816801	62.60 €			1016801	62.60 €	1216801	99.50 €
26,50 mm.#							0816601	62.60 €			1016601	62.60 €	1216601	99.50 €
25,00 mm.#					0716401	62.60 €	0816401	62.60 €	0916401	62.60 €	1016401	62.60 €	1216401	99.50 €
22,40 mm.#					0716201	62.60 €	0816201	62.60 €	0916201	62.60 €	1016201	62.60 €	1216201	99.50 €
20,00 mm.#					0716001	62.60 €	0816001	62.60 €	0916001	62.60 €	1016001	62.60 €	1216001	99.50 €
19,00 mm.#					0715801	59.50 €	0815801	59.50 €	0915801	59.50 €	1015801	59.50 €	1215801	97.00 €
18,00 mm.#					0715601	59.50 €	0815601	59.50 €	0915601	59.50 €	1015601	59.50 €	1215601	97.00 €
16,00 mm.#					0715401	59.50 €	0815401	59.50 €	0915401	59.50 €	1015401	59.50 €	1215401	97.00 €
14,00 mm.#					0715201	59.50 €	0815201	59.50 €	0915201	59.50 €	1015201	59.50 €	1215201	97.00 €
13,20 mm.#					0715001	59.50 €	0815001	59.50 €	0915001	59.50 €	1015001	59.50 €	1215001	97.00 €
12,50 mm.#					0714801	59.50 €	0814801	59.50 €	0914801	59.50 €	1014801	59.50 €	1214801	97.00 €
11,20 mm.#					0714601	59.50 €	0814601	59.50 €	0914601	59.50 €	1014601	59.50 €	1214601	97.00 €
10,00 mm.#					0714401	59.50 €	0814401	59.50 €	0914401	59.50 €	1014401	59.50 €	1214401	97.00 €

**TEST SIEVES IN STAINLESS STEEL
ACCORDING ISO 3310.1
WIRE MESH**

Abertura	200 mm.Ø Air Jet 25 mm.H		200 mm.Ø New System 25 mm.H		200 mm.Ø x 25 mm.H		200 mm.Ø x 50 mm.H		203 mm.Ø x 25 mm.H		203 mm.Ø x 50 mm.H		250 mm.Ø x 70 mm.H	
9,50 mm.#					0714201	59.50 €	0814201	59.50 €	0914201	59.50 €	1014201	59.50 €	1214201	97.00 €
9,00 mm.#					0714001	59.50 €	0814001	59.50 €	0914001	59.50 €	1014001	59.50 €	1214001	97.00 €
8,00 mm.#					0713801	59.50 €	0813801	59.50 €	0913801	59.50 €	1013801	59.50 €	1213801	97.00 €
7,10 mm.#					0713601	59.50 €	0813601	59.50 €	0913601	59.50 €	1013601	59.50 €	1213601	97.00 €
6,70 mm.#					0713401	59.50 €	0813401	59.50 €	0913401	59.50 €	1013401	59.50 €	1213401	97.00 €
6,30 mm.#					0713201	59.50 €	0813201	59.50 €	0913201	59.50 €	1013201	59.50 €	1213201	97.00 €
5,60 mm.#					0713001	59.50 €	0813001	59.50 €	0913001	59.50 €	1013001	59.50 €	1213001	97.00 €
5,00 mm.#					0712801	59.50 €	0812801	59.50 €	0912801	59.50 €	1012801	59.50 €	1212801	97.00 €
4,75 mm.#					0712601	59.50 €	0812601	59.50 €	0912601	59.50 €	1012601	59.50 €	1212601	97.00 €
4,50 mm.#					0712401	59.50 €	0812401	59.50 €	0912401	59.50 €	1012401	59.50 €	1212401	97.00 €
4,00 mm.#	0612201	80.85 €	0612299	94.50 €	0712201	57.50 €	0812201	57.50 €	0912201	57.50 €	1012201	57.50 €	1212201	94.50 €
3,55 mm.#	0612001	80.85 €	0612099	94.50 €	0712001	57.50 €	0812001	57.50 €	0912001	57.50 €	1012001	57.50 €	1212001	94.50 €
3,35 mm.#	0611801	80.85 €	0611899	94.50 €	0711801	57.50 €	0811801	57.50 €	0911801	57.50 €	1011801	57.50 €	1211801	94.50 €
3,15 mm.#	0611601	80.85 €	0611699	94.50 €	0711601	57.50 €	0811601	57.50 €	0911601	57.50 €	1011601	57.50 €	1211601	94.50 €
2,80 mm.#	0611401	80.85 €	0611499	94.50 €	0711401	57.50 €	0811401	57.50 €	0911401	57.50 €	1011401	57.50 €	1211401	94.50 €
2,50 mm.#	0611201	80.85 €	0611299	94.50 €	0711201	57.50 €	0811201	57.50 €	0911201	57.50 €	1011201	57.50 €	1211201	94.50 €
2,36 mm.#	0611001	80.85 €	0611099	94.50 €	0711001	57.50 €	0811001	57.50 €	0911001	57.50 €	1011001	57.50 €	1211001	94.50 €
2,24 mm.#	0610801	80.85 €	0610899	94.50 €	0710801	57.50 €	0810801	57.50 €	0910801	57.50 €	1010801	57.50 €	1210801	94.50 €
2,00 mm.#	0610601	80.85 €	0610699	94.50 €	0710601	57.50 €	0810601	57.50 €	0910601	57.50 €	1010601	57.50 €	1210601	94.50 €
1,80 mm.#	0610401	80.85 €	0610499	94.50 €	0710401	57.50 €	0810401	57.50 €	0910401	57.50 €	1010401	57.50 €	1210401	94.50 €
1,70 mm.#	0610201	80.85 €	0610299	94.50 €	0710201	57.50 €	0810201	57.50 €	0910201	57.50 €	1010201	57.50 €	1210201	94.50 €
1,60 mm.#	0610001	80.85 €	0610099	94.50 €	0710001	57.50 €	0810001	57.50 €	0910001	57.50 €	1010001	57.50 €	1210001	94.50 €
1,40 mm.#	0609801	80.85 €	0609899	94.50 €	0709801	57.50 €	0809801	57.50 €	0909801	57.50 €	1009801	57.50 €	1209801	94.50 €
1,25 mm.#	0609601	80.85 €	0609699	94.50 €	0709601	57.50 €	0809601	57.50 €	0909601	57.50 €	1009601	57.50 €	1209601	94.50 €
1,18 mm.#	0609401	80.85 €	0609499	94.50 €	0709401	57.50 €	0809401	57.50 €	0909401	57.50 €	1009401	57.50 €	1209401	94.50 €
1,12 mm.#	0609201	80.85 €	0609299	94.50 €	0709201	57.50 €	0809201	57.50 €	0909201	57.50 €	1009201	57.50 €	1209201	94.50 €
1,00 mm.#	0609001	80.85 €	0609099	94.50 €	0709001	57.50 €	0809001	57.50 €	0909001	57.50 €	1009001	57.50 €	1209001	94.50 €
900 µm.#	0608801	80.85 €	0608899	94.50 €	0708801	57.50 €	0808801	57.50 €	0908801	57.50 €	1008801	57.50 €	1208801	94.50 €
850 µm.#	0608601	80.85 €	0608699	94.50 €	0708601	57.50 €	0808601	57.50 €	0908601	57.50 €	1008601	57.50 €	1208601	94.50 €

**TEST SIEVES IN STAINLESS STEEL
ACCORDING ISO 3310.1
WIRE MESH**

Abertura	200 mm.Ø Air Jet 25 mm.H		200 mm.Ø New System 25 mm.H		200 mm.Ø x 25 mm.H		200 mm.Ø x 50 mm.H		203 mm.Ø x 25 mm.H		203 mm.Ø x 50 mm.H		250 mm.Ø x 70 mm.H	
800 µm.#	0608401	80.85 €	0608499	94.50 €	0708401	57.50 €	0808401	57.50 €	0908401	57.50 €	1008401	57.50 €	1208401	94.50 €
710 µm.#	0608201	80.85 €	0608299	94.50 €	0708201	57.50 €	0808201	57.50 €	0908201	57.50 €	1008201	57.50 €	1208201	94.50 €
630 µm.#	0608001	80.85 €	0608099	94.50 €	0708001	57.50 €	0808001	57.50 €	0908001	57.50 €	1008001	57.50 €	1208001	94.50 €
600 µm.#	0607801	80.85 €	0607899	94.50 €	0707801	57.50 €	0807801	57.50 €	0907801	57.50 €	1007801	57.50 €	1207801	94.50 €
560 µm.#	0607601	80.85 €	0607699	94.50 €	0707601	57.50 €	0807601	57.50 €	0907601	57.50 €	1007601	57.50 €	1207601	94.50 €
500 µm.#	0607401	80.85 €	0607499	94.50 €	0707401	57.50 €	0807401	57.50 €	0907401	57.50 €	1007401	57.50 €	1207401	94.50 €
450 µm.#	0607201	80.85 €	0607299	94.50 €	0707201	57.50 €	0807201	57.50 €	0907201	57.50 €	1007201	57.50 €	1207201	94.50 €
425 µm.#	0607001	80.85 €	0607099	94.50 €	0707001	57.50 €	0807001	57.50 €	0907001	57.50 €	1007001	57.50 €	1207001	94.50 €
400 µm.#	0606801	80.85 €	0606899	94.50 €	0706801	57.50 €	0806801	57.50 €	0906801	57.50 €	1006801	57.50 €	1206801	94.50 €
355 µm.#	0606601	80.85 €	0606699	94.50 €	0706601	57.50 €	0806601	57.50 €	0906601	57.50 €	1006601	57.50 €	1206601	94.50 €
315 µm.#	0606401	80.85 €	0606499	94.50 €	0706401	57.50 €	0806401	57.50 €	0906401	57.50 €	1006401	57.50 €	1206401	94.50 €
300 µm.#	0606201	80.85 €	0606299	94.50 €	0706201	57.50 €	0806201	57.50 €	0906201	57.50 €	1006201	57.50 €	1206201	94.50 €
280 µm.#	0606001	80.85 €	0606099	94.50 €	0706001	57.50 €	0806001	57.50 €	0906001	57.50 €	1006001	57.50 €	1206001	94.50 €
250 µm.#	0605801	80.85 €	0605899	94.50 €	0705801	57.50 €	0805801	57.50 €	0905801	57.50 €	1005801	57.50 €	1205801	94.50 €
224 µm.#	0605601	80.85 €	0605699	94.50 €	0705601	57.50 €	0805601	57.50 €	0905601	57.50 €	1005601	57.50 €	1205601	94.50 €
212 µm.#	0605401	80.85 €	0605499	94.50 €	0705401	57.50 €	0805401	57.50 €	0905401	57.50 €	1005401	57.50 €	1205401	94.50 €
200 µm.#	0605201	80.85 €	0605299	94.50 €	0705201	57.50 €	0805201	57.50 €	0905201	57.50 €	1005201	57.50 €	1205201	94.50 €
180 µm.#	0605001	80.85 €	0605099	94.50 €	0705001	57.50 €	0805001	57.50 €	0905001	57.50 €	1005001	57.50 €	1205001	94.50 €
160 µm.#	0604801	80.85 €	0604899	94.50 €	0704801	57.50 €	0804801	57.50 €	0904801	57.50 €	1004801	57.50 €	1204801	94.50 €
150 µm.#	0604601	80.85 €	0604699	94.50 €	0704601	57.50 €	0804601	57.50 €	0904601	57.50 €	1004601	57.50 €	1204601	94.50 €
140 µm.#	0604401	80.85 €	0604499	94.50 €	0704401	57.50 €	0804401	57.50 €	0904401	57.50 €	1004401	57.50 €	1204401	94.50 €
125 µm.#	0604201	80.85 €	0604299	94.50 €	0704201	57.50 €	0804201	57.50 €	0904201	57.50 €	1004201	57.50 €	1204201	94.50 €
112 µm.#	0604001	80.85 €	0604099	94.50 €	0704001	57.50 €	0804001	57.50 €	0904001	57.50 €	1004001	57.50 €	1204001	94.50 €
106 µm.#	0603801	80.85 €	0603899	94.50 €	0703801	57.50 €	0803801	57.50 €	0903801	57.50 €	1003801	57.50 €	1203801	94.50 €
100 µm.#	0603601	80.85 €	0603699	94.50 €	0703601	63.75 €	0803601	63.75 €	0903601	63.75 €	1003601	63.75 €	1203601	99.30 €
90 µm.#	0603401	80.85 €	0603499	94.50 €	0703401	63.75 €	0803401	63.75 €	0903401	63.75 €	1003401	63.75 €	1203401	99.30 €
80 µm.#	0603201	80.85 €	0603299	94.50 €	0703201	63.75 €	0803201	63.75 €	0903201	63.75 €	1003201	63.75 €	1203201	99.30 €
75 µm.#	0603001	80.85 €	0603099	94.50 €	0703001	63.75 €	0803001	63.75 €	0903001	63.75 €	1003001	63.75 €	1203001	99.30 €
71 µm.#	0602801	80.85 €	0602899	94.50 €	0702801	63.75 €	0802801	63.75 €	0902801	63.75 €	1002801	63.75 €	1202801	99.30 €

**TEST SIEVES IN STAINLESS STEEL
ACCORDING ISO 3310.1
WIRE MESH**

Abertura	200 mm.Ø Air Jet 25 mm.H		200 mm.Ø New System 25 mm.H		200 mm.Ø x 25 mm.H		200 mm.Ø x 50 mm.H		203 mm.Ø x 25 mm.H		203 mm.Ø x 50 mm.H		250 mm.Ø x 70 mm.H	
	Code	Price	Code	Price	Code	Price	Code	Price	Code	Price	Code	Price	Code	Price
63 µm.#	0602601	80.85 €	0602699	94.50 €	0702601	63.75 €	0802601	63.75 €	0902601	63.75 €	1002601	63.75 €	1202601	99.30 €
56 µm.#	0602401	86.10 €	0602499	99.75 €	0702401	77.30 €	0802401	77.30 €	0902401	77.30 €	1002401	77.30 €	1202401	120.00 €
53 µm.#	0602201	86.10 €	0602299	99.75 €	0702201	77.30 €	0802201	77.30 €	0902201	77.30 €	1002201	77.30 €	1202201	120.00 €
50 µm.#	0602001	86.10 €	0602099	99.75 €	0702001	77.30 €	0802001	77.30 €	0902001	77.30 €	1002001	77.30 €	1202001	120.00 €
45 µm.#	0601801	86.10 €	0601899	99.75 €	0701801	77.30 €	0801801	77.30 €	0901801	77.30 €	1001801	77.30 €	1201801	120.00 €
40 µm.#	0601601	134.40 €	0601699	141.75 €	0701601	88.80 €	0801601	88.80 €	0901601	88.80 €	1001601	88.80 €	1201601	130.60 €
38 µm.#	0601401	134.40 €	0601499	141.75 €	0701401	88.80 €	0801401	88.80 €	0901401	88.80 €	1001401	88.80 €	1201401	130.60 €
36 µm.#	0601201	134.40 €	0601299	141.75 €	0701201	88.80 €	0801201	88.80 €	0901201	88.80 €	1001201	88.80 €	1201201	130.60 €
32 µm.#	0601001	143.85 €	0601099	157.50 €	0701001	101.50 €	0801001	101.50 €	0901001	101.50 €	1001001	101.50 €	1201001	146.30 €
25 µm.#	0600801	156.00 €	0600899	168.00 €	0700801	136.50 €	0800801	136.50 €	0900801	136.50 €	1000801	136.50 €	1200801	212.00 €
20 µm.#	0600601	410.00 €	0600699	425.00 €	0700601	370.00 €	0800601	370.00 €	0900601	370.00 €	1000601	370.00 €	1200601	

SPECIAL OPENINGS (NON STANDARDIZED)

Abertura	200 mm.Ø Air Jet 25 mm.H		200 mm.Ø New System 25 mm.H		200 mm.Ø x 25 mm.H		200 mm.Ø x 50 mm.H		203 mm.Ø x 25 mm.H		203 mm.Ø x 50 mm.H		250 mm.Ø x 70 mm.H	
	Code	Price	Code	Price	Code	Price	Code	Price	Code	Price	Code	Price	Code	Price
6,00 mm.#					0713101	59.50 €	0813101	59.50 €	0913101	59.50 €	1013101	59.50 €	1213101	97.00 €
3,00 mm.#	0611501	80.85 €	0611599	94.50 €	0711501	57.50 €	0811501	57.50 €	0911501	57.50 €	1011501	57.50 €	1211501	94.50 €
1,50 mm.#	0609901	80.85 €	0609999	94.50 €	0709901	57.50 €	0809901	57.50 €	0909901	57.50 €	1009901	57.50 €	1209901	94.50 €
750 µm.#	0608301	80.85 €	0608399	94.50 €	0708301	57.50 €	0808301	57.50 €	0908301	57.50 €	1008301	57.50 €	1208301	94.50 €

**TEST SIEVES IN STAINLESS STEEL
ACCORDING ISO 3310.1
WIRE MESH**

	300 mm.Ø x 70 mm.H		305 mm.Ø x 80 mm.H		315 mm.Ø x 80 mm.H		350 mm.Ø x 60 mm.H		400 mm.Ø x 65 mm.H		400 mm.Ø x 100 mm.H		450 mm.Ø x 100 mm.H	
Opening	M-LB code	Price	M-LB code	Price	M-LB code	Price	M-LB code	Price	M-LB code	Price	M-LB code	Price	M-LB code	Price
125,00 mm.#	1620201	131.00 €	1720201	131.00 €	1520201	131.00 €	3220201	161.00 €	2420201	207.00 €	1820201	217.00 €	1920201	280.00 €
112,00 mm.#	1620001	131.00 €	1720001	131.00 €	1520001	131.00 €	3220001	161.00 €	2420001	207.00 €	1820001	217.00 €	1920001	280.00 €
106,00 mm.#	1619801	131.00 €	1719801	131.00 €	1519801	131.00 €	3219801	161.00 €	2419801	207.00 €	1819801	217.00 €	1919801	280.00 €
100,00 mm.#	1619601	131.00 €	1719601	131.00 €	1519601	131.00 €	3219601	161.00 €	2419601	207.00 €	1819601	217.00 €	1919601	280.00 €
90,00 mm.#	1619401	131.00 €	1719401	131.00 €	1519401	131.00 €	3219401	161.00 €	2419401	207.00 €	1819401	217.00 €	1919401	280.00 €
80,00 mm.#	1619201	131.00 €	1719201	131.00 €	1519201	131.00 €	3219201	161.00 €	2419201	207.00 €	1819201	217.00 €	1919201	280.00 €
75,00 mm.#	1619001	131.00 €	1719001	131.00 €	1519001	131.00 €	3219001	161.00 €	2419001	207.00 €	1819001	217.00 €	1919001	280.00 €
71,00 mm.#	1618801	131.00 €	1718801	131.00 €	1518801	131.00 €	3218801	161.00 €	2418801	207.00 €	1818801	217.00 €	1918801	280.00 €
63,00 mm.#	1618601	131.00 €	1718601	131.00 €	1518601	131.00 €	3218601	161.00 €	2418601	207.00 €	1818601	217.00 €	1918601	280.00 €
56,00 mm.#	1618401	131.00 €	1718401	131.00 €	1518401	131.00 €	3218401	161.00 €	2418401	207.00 €	1818401	217.00 €	1918401	280.00 €
53,00 mm.#	1618201	131.00 €	1718201	131.00 €	1518201	131.00 €	3218201	161.00 €	2418201	207.00 €	1818201	217.00 €	1918201	280.00 €
50,00 mm.#	1618001	131.00 €	1718001	131.00 €	1518001	131.00 €	3218001	161.00 €	2418001	207.00 €	1818001	217.00 €	1918001	280.00 €
45,00 mm.#	1617801	131.00 €	1717801	131.00 €	1517801	131.00 €	3217801	161.00 €	2417801	207.00 €	1817801	217.00 €	1917801	280.00 €
40,00 mm.#	1617601	131.00 €	1717601	131.00 €	1517601	131.00 €	3217601	161.00 €	2417601	207.00 €	1817601	217.00 €	1917601	280.00 €
37,50 mm.#	1617401	131.00 €	1717401	131.00 €	1517401	131.00 €	3217401	161.00 €	2417401	207.00 €	1817401	217.00 €	1917401	280.00 €
35,50mm.#	1617201	131.00 €	1717201	131.00 €	1517201	131.00 €	3217201	161.00 €	2417201	207.00 €	1817201	217.00 €	1917201	280.00 €
31,50 mm.#	1617001	131.00 €	1717001	131.00 €	1517001	131.00 €	3217001	161.00 €	2417001	207.00 €	1817001	217.00 €	1917001	280.00 €
28,00 mm.#	1616801	131.00 €	1716801	131.00 €	1516801	131.00 €	3216801	161.00 €	2416801	207.00 €	1816801	217.00 €	1916801	280.00 €
26,50 mm.#	1616601	131.00 €	1716601	131.00 €	1516601	131.00 €	3216601	161.00 €	2416601	207.00 €	1816601	217.00 €	1916601	280.00 €
25,00 mm.#	1616401	131.00 €	1716401	131.00 €	1516401	131.00 €	3216401	161.00 €	2416401	207.00 €	1816401	217.00 €	1916401	280.00 €
22,40 mm.#	1616201	131.00 €	1716201	131.00 €	1516201	131.00 €	3216201	161.00 €	2416201	207.00 €	1816201	217.00 €	1916201	280.00 €
20,00 mm.#	1616001	131.00 €	1716001	131.00 €	1516001	131.00 €	3216001	161.00 €	2416001	207.00 €	1816001	217.00 €	1916001	280.00 €
19,00 mm.#	1615801	114.90 €	1715801	114.90 €	1515801	114.90 €	3215801	145.00 €	2415801	182.00 €	1815801	192.00 €	1915801	254.00 €
18,00 mm.#	1615601	114.90 €	1715601	114.90 €	1515601	114.90 €	3215601	145.00 €	2415601	182.00 €	1815601	192.00 €	1915601	254.00 €
16,00 mm.#	1615401	114.90 €	1715401	114.90 €	1515401	114.90 €	3215401	145.00 €	2415401	182.00 €	1815401	192.00 €	1915401	254.00 €
14,00 mm.#	1615201	114.90 €	1715201	114.90 €	1515201	114.90 €	3215201	145.00 €	2415201	182.00 €	1815201	192.00 €	1915201	254.00 €
13,20 mm.#	1615001	114.90 €	1715001	114.90 €	1515001	114.90 €	3215001	145.00 €	2415001	182.00 €	1815001	192.00 €	1915001	254.00 €
12,50 mm.#	1614801	114.90 €	1714801	114.90 €	1514801	114.90 €	3214801	145.00 €	2414801	182.00 €	1814801	192.00 €	1914801	254.00 €
11,20 mm.#	1614601	114.90 €	1714601	114.90 €	1514601	114.90 €	3214601	145.00 €	2414601	182.00 €	1814601	192.00 €	1914601	254.00 €
10,00 mm.#	1614401	114.90 €	1714401	114.90 €	1514401	114.90 €	3214401	145.00 €	2414401	182.00 €	1814401	192.00 €	1914401	254.00 €

**TEST SIEVES IN STAINLESS STEEL
ACCORDING ISO 3310.1
WIRE MESH**

Abertura	300 mm.Ø x 70 mm.H		305 mm.Ø x 80 mm.H		315 mm.Ø x 80 mm.H		350 mm.Ø x 70 mm.H		400 mm.Ø x 65 mm.H		400 mm.Ø x 100 mm.H		450 mm.Ø x 100 mm.H	
9,50 mm.#	1614201	114.90 €	1714201	114.90 €	1514201	114.90 €	3214201	145.00 €	2414201	182.00 €	1814201	192.00 €	1914201	254.00 €
9,00 mm.#	1614001	114.90 €	1714001	114.90 €	1514001	114.90 €	3214001	145.00 €	2414001	182.00 €	1814001	192.00 €	1914001	254.00 €
8,00 mm.#	1613801	114.90 €	1713801	114.90 €	1513801	114.90 €	3213801	145.00 €	2413801	182.00 €	1813801	192.00 €	1913801	254.00 €
7,10 mm.#	1613601	114.90 €	1713601	114.90 €	1513601	114.90 €	3213601	145.00 €	2413601	182.00 €	1813601	192.00 €	1913601	254.00 €
6,70 mm.#	1613401	114.90 €	1713401	114.90 €	1513401	114.90 €	3213401	145.00 €	2413401	182.00 €	1813401	192.00 €	1913401	254.00 €
6,30 mm.#	1613201	114.90 €	1713201	114.90 €	1513201	114.90 €	3213201	145.00 €	2413201	182.00 €	1813201	192.00 €	1913201	254.00 €
5,60 mm.#	1613001	114.90 €	1713001	114.90 €	1513001	114.90 €	3213001	145.00 €	2413001	182.00 €	1813001	192.00 €	1913001	254.00 €
5,00 mm.#	1612801	114.90 €	1712801	114.90 €	1512801	114.90 €	3212801	145.00 €	2412801	182.00 €	1812801	192.00 €	1912801	254.00 €
4,75 mm.#	1612601	114.90 €	1712601	114.90 €	1512601	114.90 €	3212601	145.00 €	2412601	182.00 €	1812601	192.00 €	1912601	254.00 €
4,50 mm.#	1612401	114.90 €	1712401	114.90 €	1512401	114.90 €	3212401	145.00 €	2412401	182.00 €	1812401	192.00 €	1912401	254.00 €
4,00 mm.#	1612201	99.30 €	1712201	99.30 €	1512201	99.30 €	3212201	129.00 €	2412201	175.00 €	1812201	185.00 €	1912201	245.00 €
3,55 mm.#	1612001	99.30 €	1712001	99.30 €	1512001	99.30 €	3212001	129.00 €	2412001	175.00 €	1812001	185.00 €	1912001	245.00 €
3,35 mm.#	1611801	99.30 €	1711801	99.30 €	1511801	99.30 €	3211801	129.00 €	2411801	175.00 €	1811801	185.00 €	1911801	245.00 €
3,15 mm.#	1611601	99.30 €	1711601	99.30 €	1511601	99.30 €	3211601	129.00 €	2411601	175.00 €	1811601	185.00 €	1911601	245.00 €
2,80 mm.#	1611401	99.30 €	1711401	99.30 €	1511401	99.30 €	3211401	129.00 €	2411401	175.00 €	1811401	185.00 €	1911401	245.00 €
2,50 mm.#	1611201	99.30 €	1711201	99.30 €	1511201	99.30 €	3211201	129.00 €	2411201	175.00 €	1811201	185.00 €	1911201	245.00 €
2,36 mm.#	1611001	99.30 €	1711001	99.30 €	1511001	99.30 €	3211001	129.00 €	2411001	175.00 €	1811001	185.00 €	1911001	245.00 €
2,24 mm.#	1610801	99.30 €	1710801	99.30 €	1510801	99.30 €	3210801	129.00 €	2410801	175.00 €	1810801	185.00 €	1910801	245.00 €
2,00 mm.#	1610601	99.30 €	1710601	99.30 €	1510601	99.30 €	3210601	129.00 €	2410601	175.00 €	1810601	185.00 €	1910601	245.00 €
1,80 mm.#	1610401	99.30 €	1710401	99.30 €	1510401	99.30 €	3210401	129.00 €	2410401	175.00 €	1810401	185.00 €	1910401	245.00 €
1,70 mm.#	1610201	99.30 €	1710201	99.30 €	1510201	99.30 €	3210201	129.00 €	2410201	175.00 €	1810201	185.00 €	1910201	245.00 €
1,60 mm.#	1610001	99.30 €	1710001	99.30 €	1510001	99.30 €	3210001	129.00 €	2410001	175.00 €	1810001	185.00 €	1910001	245.00 €
1,40 mm.#	1609801	99.30 €	1709801	99.30 €	1509801	99.30 €	3209801	129.00 €	2409801	175.00 €	1809801	185.00 €	1909801	245.00 €
1,25 mm.#	1609601	99.30 €	1709601	99.30 €	1509601	99.30 €	3209601	129.00 €	2409601	175.00 €	1809601	185.00 €	1909601	245.00 €
1,18 mm.#	1609401	99.30 €	1709401	99.30 €	1509401	99.30 €	3209401	129.00 €	2409401	175.00 €	1809401	185.00 €	1909401	245.00 €
1,12 mm.#	1609201	99.30 €	1709201	99.30 €	1509201	99.30 €	3209201	129.00 €	2409201	175.00 €	1809201	185.00 €	1909201	245.00 €
1,00 mm.#	1609001	99.30 €	1709001	99.30 €	1509001	99.30 €	3209001	129.00 €	2409001	175.00 €	1809001	185.00 €	1909001	245.00 €
900 µm.#	1608801	99.30 €	1708801	99.30 €	1508801	99.30 €	3208801	129.00 €	2408801	175.00 €	1808801	185.00 €	1908801	245.00 €
850 µm.#	1608601	99.30 €	1708601	99.30 €	1508601	99.30 €	3208601	129.00 €	2408601	175.00 €	1808601	185.00 €	1908601	245.00 €

**TEST SIEVES IN STAINLESS STEEL
ACCORDING ISO 3310.1
WIRE MESH**

Abertura	300 mm.Ø x 70 mm.H		305 mm.Ø x 80 mm.H		315 mm.Ø x 80 mm.H		350 mm.Ø x 70 mm.H		400 mm.Ø x 65 mm.H		400 mm.Ø x 100 mm.H		450 mm.Ø x 100 mm.H	
	Code	Price	Code	Price	Code	Price	Code	Price	Code	Price	Code	Price	Code	Price
800 µm.#	1608401	99.30 €	1708401	99.30 €	1508401	99.30 €	3208401	129.00 €	2408401	175.00 €	1808401	185.00 €	1908401	245.00 €
710 µm.#	1608201	99.30 €	1708201	99.30 €	1508201	99.30 €	3208201	129.00 €	2408201	175.00 €	1808201	185.00 €	1908201	245.00 €
630 µm.#	1608001	99.30 €	1708001	99.30 €	1508001	99.30 €	3208001	129.00 €	2408001	175.00 €	1808001	185.00 €	1908001	245.00 €
600 µm.#	1607801	99.30 €	1707801	99.30 €	1507801	99.30 €	3207801	129.00 €	2407801	175.00 €	1807801	185.00 €	1907801	245.00 €
560 µm.#	1607601	99.30 €	1707601	99.30 €	1507601	99.30 €	3207601	129.00 €	2407601	175.00 €	1807601	185.00 €	1907601	245.00 €
500 µm.#	1607401	99.30 €	1707401	99.30 €	1507401	99.30 €	3207401	129.00 €	2407401	175.00 €	1807401	185.00 €	1907401	245.00 €
450 µm.#	1607201	99.30 €	1707201	99.30 €	1507201	99.30 €	3207201	129.00 €	2407201	175.00 €	1807201	185.00 €	1907201	245.00 €
425 µm.#	1607001	99.30 €	1707001	99.30 €	1507001	99.30 €	3207001	129.00 €	2407001	175.00 €	1807001	185.00 €	1907001	245.00 €
400 µm.#	1606801	99.30 €	1706801	99.30 €	1506801	99.30 €	3206801	129.00 €	2406801	175.00 €	1806801	185.00 €	1906801	245.00 €
355 µm.#	1606601	99.30 €	1706601	99.30 €	1506601	99.30 €	3206601	129.00 €	2406601	175.00 €	1806601	185.00 €	1906601	245.00 €
315 µm.#	1606401	99.30 €	1706401	99.30 €	1506401	99.30 €	3206401	129.00 €	2406401	175.00 €	1806401	185.00 €	1906401	245.00 €
300 µm.#	1606201	99.30 €	1706201	99.30 €	1506201	99.30 €	3206201	129.00 €	2406201	175.00 €	1806201	185.00 €	1906201	245.00 €
280 µm.#	1606001	99.30 €	1706001	99.30 €	1506001	99.30 €	3206001	129.00 €	2406001	175.00 €	1806001	185.00 €	1906001	245.00 €
250 µm.#	1605801	99.30 €	1705801	99.30 €	1505801	99.30 €	3205801	129.00 €	2405801	175.00 €	1805801	185.00 €	1905801	245.00 €
224 µm.#	1605601	99.30 €	1705601	99.30 €	1505601	99.30 €	3205601	129.00 €	2405601	175.00 €	1805601	185.00 €	1905601	245.00 €
212 µm.#	1605401	99.30 €	1705401	99.30 €	1505401	99.30 €	3205401	129.00 €	2405401	175.00 €	1805401	185.00 €	1905401	245.00 €
200 µm.#	1605201	99.30 €	1705201	99.30 €	1505201	99.30 €	3205201	129.00 €	2405201	175.00 €	1805201	185.00 €	1905201	245.00 €
180 µm.#	1605001	99.30 €	1705001	99.30 €	1505001	99.30 €	3205001	129.00 €	2405001	175.00 €	1805001	185.00 €	1905001	245.00 €
160 µm.#	1604801	99.30 €	1704801	99.30 €	1504801	99.30 €	3204801	129.00 €	2404801	175.00 €	1804801	185.00 €	1904801	245.00 €
150 µm.#	1604601	99.30 €	1704601	99.30 €	1504601	99.30 €	3204601	129.00 €	2404601	175.00 €	1804601	185.00 €	1904601	245.00 €
140 µm.#	1604401	99.30 €	1704401	99.30 €	1504401	99.30 €	3204401	129.00 €	2404401	175.00 €	1804401	185.00 €	1904401	245.00 €
125 µm.#	1604201	99.30 €	1704201	99.30 €	1504201	99.30 €	3204201	129.00 €	2404201	175.00 €	1804201	185.00 €	1904201	245.00 €
112 µm.#	1604001	99.30 €	1704001	99.30 €	1504001	99.30 €	3204001	129.00 €	2404001	175.00 €	1804001	185.00 €	1904001	245.00 €
106 µm.#	1603801	99.30 €	1703801	99.30 €	1503801	99.30 €	3203801	129.00 €	2403801	175.00 €	1803801	185.00 €	1903801	245.00 €
100 µm.#	1603601	110.00 €	1703601	110.00 €	1503601	110.00 €	3203601	142.00 €	2403601	181.00 €	1803601	191.00 €	1903601	254.00 €
90 µm.#	1603401	110.00 €	1703401	110.00 €	1503401	110.00 €	3203401	142.00 €	2403401	181.00 €	1803401	191.00 €	1903401	254.00 €
80 µm.#	1603201	110.00 €	1703201	110.00 €	1503201	110.00 €	3203201	142.00 €	2403201	181.00 €	1803201	191.00 €	1903201	254.00 €
75 µm.#	1603001	110.00 €	1703001	110.00 €	1503001	110.00 €	3203001	142.00 €	2403001	181.00 €	1803001	191.00 €	1903001	254.00 €
71 µm.#	1602801	110.00 €	1702801	110.00 €	1502801	110.00 €	3202801	142.00 €	2402801	181.00 €	1802801	191.00 €	1902801	254.00 €

**TEST SIEVES IN STAINLESS STEEL
ACCORDING ISO 3310.1
WIRE MESH**

Abertura	300 mm.Ø x 70 mm.H		305 mm.Ø x 80 mm.H		315 mm.Ø x 80 mm.H		350 mm.Ø x 70 mm.H		400 mm.Ø x 65 mm.H		400 mm.Ø x 100 mm.H		450 mm.Ø x 100 mm.H	
	Code	Price	Code	Price	Code	Price	Code	Price	Code	Price	Code	Price	Code	Price
63 µm.#	1602601	110.00 €	1702601	110.00 €	1502601	110.00 €	3202601	142.00 €	2402601	181.00 €	1802601	191.00 €	1902601	254.00 €
56 µm.#	1602401	130.00 €	1702401	130.00 €	1502401	130.00 €	3202401	165.00 €	2402401	220.00 €	1802401	230.00 €	1902401	310.00 €
53 µm.#	1602201	130.00 €	1702201	130.00 €	1502201	130.00 €	3202201	165.00 €	2402201	220.00 €	1802201	230.00 €	1902201	310.00 €
50 µm.#	1602001	130.00 €	1702001	130.00 €	1502001	130.00 €	3202001	165.00 €	2402001	220.00 €	1802001	230.00 €	1902001	310.00 €
45 µm.#	1601801	130.00 €	1701801	130.00 €	1501801	130.00 €	3201801	165.00 €	2401801	220.00 €	1801801	230.00 €	1901801	310.00 €
40 µm.#	1601601	175.00 €	1701601	175.00 €	1501601	175.00 €	3201601	220.00 €	2401601	270.00 €	1801601	280.00 €	1901601	380.00 €
38 µm.#	1601401	175.00 €	1701401	175.00 €	1501401	175.00 €	3201401	220.00 €	2401401	270.00 €	1801401	280.00 €	1901401	380.00 €
36 µm.#	1601201	175.00 €	1701201	175.00 €	1501201	175.00 €	3201201	220.00 €	2401201	270.00 €	1801201	280.00 €	1901201	380.00 €
32 µm.#	1601001	180.00 €	1701001	180.00 €	1501001	180.00 €	3201001	225.00 €	2401001	300.00 €	1801001	310.00 €	1901001	420.00 €
25 µm.#	1600801	255.00 €	1700801	255.00 €	1500801	255.00 €	3200801	275.00 €	2400801	400.00 €	1800801	410.00 €	1900801	480.00 €
20 µm.#	1600601		1700601		1500601		3200601		2400601		1800601		1900601	

SPECIAL OPENINGS (NON STANDARDIZED)

Abertura	300 mm.Ø x 70 mm.H		305 mm.Ø x 80 mm.H		315 mm.Ø x 80 mm.H		350 mm.Ø x 70 mm.H		400 mm.Ø x 65 mm.H		400 mm.Ø x 100 mm.H		450 mm.Ø x 100 mm.H	
	Code	Price	Code	Price	Code	Price	Code	Price	Code	Price	Code	Price	Code	Price
6,00 mm.#	1613101	114.90 €	1713101	114.90 €	1513101	114.90 €	3213101	145.00 €	2413101	182.00 €	1813101	192.00 €	1913101	254.00 €
3,00 mm.#	1611501	99.30 €	1711501	99.30 €	1511501	99.30 €	3211501	129.00 €	2411501	175.00 €	1811501	185.00 €	1911501	245.00 €
1,50 mm.#	1609901	99.30 €	1709901	99.30 €	1509901	99.30 €	3209901	129.00 €	2409901	175.00 €	1809901	185.00 €	1909901	245.00 €
750 µm.#	1608301	99.30 €	1708301	99.30 €	1508301	99.30 €	3208301	129.00 €	2408301	175.00 €	1808301	185.00 €	1908301	245.00 €

**STAINLESS STEEL TEST SIEVES
ACCORDING ISO 3310.2
PERFORATED PLATES - SQUARE OPENINGS**

Opening	200 mm.Ø x 50 mm.H		203 mm.Ø x 50 mm.H		300 mm.Ø x 70 mm.H		305 mm.Ø x 80 mm.H		350 mm.Ø x 60 mm.H		400 mm.Ø x 100 mm.H		450 mm.Ø x 100 mm.H	
	M-LB code	Price	M-LB code	Price	M-LB code	Price	M-LB code	Price	M-LB code	Price	M-LB code	Price	M-LB code	Price
125,00 mm.■	0820261	75.75 €	1020261	75.75 €	1620261	113.90 €	1720261	113.90 €	3220261	170.00 €	1820261	192.00 €	1920261	240.00 €
112,00 mm.■	0820061	75.75 €	1020061	75.75 €	1620061	113.90 €	1720061	113.90 €	3220061	170.00 €	1820061	192.00 €	1920061	240.00 €
106,00 mm.■	0819861	75.75 €	1019861	75.75 €	1619861	113.90 €	1719861	113.90 €	3219861	170.00 €	1819861	192.00 €	1919861	240.00 €
100,00 mm.■	0819661	75.75 €	1019661	75.75 €	1619661	113.90 €	1719661	113.90 €	3219661	170.00 €	1819661	192.00 €	1919661	240.00 €
90,00 mm.■	0819461	75.75 €	1019461	75.75 €	1619461	113.90 €	1719461	113.90 €	3219461	170.00 €	1819461	192.00 €	1919461	240.00 €
80,00 mm.■	0819261	75.75 €	1019261	75.75 €	1619261	113.90 €	1719261	113.90 €	3219261	170.00 €	1819261	192.00 €	1919261	240.00 €
75,00 mm.■	0819061	75.75 €	1019061	75.75 €	1619061	113.90 €	1719061	113.90 €	3219061	170.00 €	1819061	192.00 €	1919061	240.00 €
71,00 mm.■	0818861	75.75 €	1018861	75.75 €	1618861	113.90 €	1718861	113.90 €	3218861	170.00 €	1818861	192.00 €	1918861	240.00 €
63,00 mm.■	0818661	75.75 €	1018661	75.75 €	1618661	113.90 €	1718661	113.90 €	3218661	170.00 €	1818661	192.00 €	1918661	240.00 €
56,00 mm.■	0818461	75.75 €	1018461	75.75 €	1618461	113.90 €	1718461	113.90 €	3218461	170.00 €	1818461	192.00 €	1918461	240.00 €
53,00 mm.■	0818261	75.75 €	1018261	75.75 €	1618261	113.90 €	1718261	113.90 €	3218261	170.00 €	1818261	192.00 €	1918261	240.00 €
50,00 mm.■	0818061	75.75 €	1018061	75.75 €	1618061	113.90 €	1718061	113.90 €	3218061	170.00 €	1818061	192.00 €	1918061	240.00 €
45,00 mm.■	0817861	75.75 €	1017861	75.75 €	1617861	113.90 €	1717861	113.90 €	3217861	170.00 €	1817861	192.00 €	1917861	240.00 €
40,00 mm.■	0817661	75.75 €	1017661	75.75 €	1617661	113.90 €	1717661	113.90 €	3217661	170.00 €	1817661	192.00 €	1917661	240.00 €
37,50 mm.■	0817461	75.75 €	1017461	75.75 €	1617461	113.90 €	1717461	113.90 €	3217461	170.00 €	1817461	192.00 €	1917461	240.00 €
35,50 mm.■	0817261	75.75 €	1017261	75.75 €	1617261	113.90 €	1717261	113.90 €	3217261	170.00 €	1817261	192.00 €	1917261	240.00 €
31,50 mm.■	0817061	75.75 €	1017061	75.75 €	1617061	113.90 €	1717061	113.90 €	3217061	170.00 €	1817061	192.00 €	1917061	240.00 €
28,00 mm.■	0816861	75.75 €	1016861	75.75 €	1616861	113.90 €	1716861	113.90 €	3216861	170.00 €	1816861	192.00 €	1916861	240.00 €
26,50 mm.■	0816661	75.75 €	1016661	75.75 €	1616661	113.90 €	1716661	113.90 €	3216661	170.00 €	1816661	192.00 €	1916661	240.00 €
25,00 mm.■	0816461	75.75 €	1016461	75.75 €	1616461	113.90 €	1716461	113.90 €	3216461	170.00 €	1816461	192.00 €	1916461	240.00 €
22,40 mm.■	0816261	75.75 €	1016261	75.75 €	1616261	113.90 €	1716261	113.90 €	3216261	170.00 €	1816261	192.00 €	1916261	240.00 €
20,00 mm.■	0816061	75.75 €	1016061	75.75 €	1616061	113.90 €	1716061	113.90 €	3216061	170.00 €	1816061	192.00 €	1916061	240.00 €
19,00 mm.■	0815861	75.75 €	1015861	75.75 €	1615861	113.90 €	1715861	113.90 €	3215861	170.00 €	1815861	192.00 €	1915861	240.00 €
18,00 mm.■	0815661	75.75 €	1015661	75.75 €	1615661	113.90 €	1715661	113.90 €	3215661	170.00 €	1815661	192.00 €	1915661	240.00 €
16,00 mm.■	0815461	79.95 €	1015461	79.95 €	1615461	129.50 €	1715461	129.50 €	3215461	156.00 €	1815461	177.65 €	1915461	225.00 €
14,00 mm.■	0815261	79.95 €	1015261	79.95 €	1615261	129.50 €	1715261	129.50 €	3215261	156.00 €	1815261	177.65 €	1915261	225.00 €
13,20 mm.■	0815061	79.95 €	1015061	79.95 €	1615061	129.50 €	1715061	129.50 €	3215061	156.00 €	1815061	177.65 €	1915061	225.00 €
12,50 mm.■	0814861	79.95 €	1014861	79.95 €	1614861	129.50 €	1714861	129.50 €	3214861	156.00 €	1814861	177.65 €	1914861	225.00 €
11,20 mm.■	0814661	79.95 €	1014661	79.95 €	1614661	129.50 €	1714661	129.50 €	3214661	156.00 €	1814661	177.65 €	1914661	225.00 €

**STAINLESS STEEL TEST SIEVES
ACCORDING ISO 3310.2
PERFORATED PLATES - SQUARE OPENINGS**

Opening	200 mm.Ø x 50 mm.H		203 mm.Ø x 50 mm.H		300 mm.Ø x 70 mm.H		305 mm.Ø x 80 mm.H		350 mm.Ø x 60 mm.H		400 mm.Ø x 100 mm.H		450 mm.Ø x 100 mm.H	
	M-LB code	Price	M-LB code	Price	M-LB code	Price	M-LB code	Price	M-LB code	Price	M-LB code	Price	M-LB code	Price
10,00 mm. ■	0814461	79.95 €	1014461	79.95 €	1614461	129.50 €	1714461	129.50 €	3214461	156.00 €	1814461	177.65 €	1914461	225.00 €
9,50 mm. ■	0814261	79.95 €	1014261	79.95 €	1614261	129.50 €	1714261	129.50 €	3214261	156.00 €	1814261	177.65 €	1914261	225.00 €
9,00 mm. ■	0814061	79.95 €	1014061	79.95 €	1614061	129.50 €	1714061	129.50 €	3214061	156.00 €	1814061	177.65 €	1914061	225.00 €
8,00 mm. ■	0813861	79.95 €	1013861	79.95 €	1613861	129.50 €	1713861	129.50 €	3213861	156.00 €	1813861	177.65 €	1913861	225.00 €
7,10 mm. ■	0813661	79.95 €	1013661	79.95 €	1613661	129.50 €	1713661	129.50 €	3213661	156.00 €	1813661	177.65 €	1913661	225.00 €
6,70 mm. ■	0813461	79.95 €	1013461	79.95 €	1613461	129.50 €	1713461	129.50 €	3213461	156.00 €	1813461	177.65 €	1913461	225.00 €
6,30 mm. ■	0813261	79.95 €	1013261	79.95 €	1613261	129.50 €	1713261	129.50 €	3213261	180.00 €	1813261	209.00 €	1913261	282.00 €
5,60 mm. ■	0813061	79.95 €	1013061	79.95 €	1613061	129.50 €	1713061	129.50 €	3213061	180.00 €	1813061	209.00 €	1913061	282.00 €
5,00 mm. ■	0812861	79.95 €	1012861	79.95 €	1612861	129.50 €	1712861	129.50 €	3212861	180.00 €	1812861	209.00 €	1912861	282.00 €
4,75 mm. ■	0812661	79.95 €	1012661	79.95 €	1612661	129.50 €	1712661	129.50 €	3212661	180.00 €	1812661	209.00 €	1912661	282.00 €
4,50 mm. ■	0812461	79.95 €	1012461	79.95 €	1612461	129.50 €	1712461	129.50 €	3212461	180.00 €	1812461	209.00 €	1912461	282.00 €
4,00 mm. ■	0812261	79.95 €	1012261	79.95 €	1612261	129.50 €	1712261	129.50 €	3212261	180.00 €	1812261	209.00 €	1912261	282.00 €

NON STANDARDIZED OPENINGS

Opening	M-LB code	Price	M-LB code	Price	M-LB code	Price	M-LB code	Price	M-LB code	Price	M-LB code	Price	M-LB code	Price
150,00 mm. ■	0820561	75.75 €	1020561	75.75 €	1620561	113.90 €	1720561	113.90 €	3220561	170.00 €	1820561	192.00 €	1920561	240.00 €
55,00 mm. ■	0818361	75.75 €	1018361	75.75 €	1618361	113.90 €	1718361	113.90 €	3218361	170.00 €	1818361	192.00 €	1918361	240.00 €
48,00 mm. ■	0817961	75.75 €	1017961	75.75 €	1617961	113.90 €	1717961	113.90 €	3217961	170.00 €	1817961	192.00 €	1917961	240.00 €
35,00 mm. ■	0817181	75.75 €	1017181	75.75 €	1617181	113.90 €	1717181	113.90 €	3217181	170.00 €	1817181	177.65 €	1917181	225.00 €
32,00 mm. ■	0817161	75.75 €	1017161	75.75 €	1617161	113.90 €	1717161	113.90 €	3217161	170.00 €	1817161	177.65 €	1917161	225.00 €
31,00 mm. ■	0816961	75.75 €	1016961	75.75 €	1616961	113.90 €	1716961	113.90 €	3216961	170.00 €	1816961	177.65 €	1916961	225.00 €
15,00 mm. ■	0815361	79.95 €	1015361	79.95 €	1615361	129.50 €	1715361	129.50 €	3215361	156.00 €	1815361	177.65 €	1915361	225.00 €
12,00 mm. ■	0814761	79.95 €	1014761	79.95 €	1614761	129.50 €	1714761	129.50 €	3214761	156.00 €	1814761	177.65 €	1914761	225.00 €
6,00 mm. ■	0813161	79.95 €	1013161	79.95 €	1613161	129.50 €	1713161	129.50 €	3213161	180.00 €	1813161	209.00 €	1913161	282.00 €

**STAINLESS STEEL TEST SIEVES
ACCORDING ISO 3310.2
PERFORATED PLATES - ROUND OPENINGS**

Opening	200 mm.Ø x 50 mm.H		203 mm.Ø x 50 mm.H		300 mm.Ø x 70 mm.H		305 mm.Ø x 80 mm.H		350 mm.Ø x 60 mm.H		400 mm.Ø x 100 mm.H		450 mm.Ø x 100 mm.H	
	M-LB code	Price	M-LB code	Price	M-LB code	Price	M-LB code	Price	M-LB code	Price	M-LB code	Price	M-LB code	Price
125,00 mm.●	0820241	75.25 €	1020241	75.25 €	1620241	109.70 €	1720241	109.70 €	3220241	170.00 €	1820241	188.00 €	1920241	240.00 €
112,00 mm.●	0820041	75.25 €	1020041	75.25 €	1620041	109.70 €	1720041	109.70 €	3220041	170.00 €	1820041	188.00 €	1920041	240.00 €
106,00 mm.●	0819841	75.25 €	1019841	75.25 €	1419841	109.70 €	1719841	109.70 €	3219841	170.00 €	1819841	188.00 €	1919841	240.00 €
100,00 mm.●	0819641	75.25 €	1019641	75.25 €	1419641	109.70 €	1719641	109.70 €	3219641	170.00 €	1819641	188.00 €	1919641	240.00 €
90,00 mm.●	0819441	75.25 €	1019441	75.25 €	1419441	109.70 €	1719441	109.70 €	3219441	170.00 €	1819441	188.00 €	1919441	240.00 €
80,00 mm.●	0819241	75.25 €	1019241	75.25 €	1419241	109.70 €	1719241	109.70 €	3219241	170.00 €	1819241	188.00 €	1919241	240.00 €
75,00 mm.●	0819041	75.25 €	1019041	75.25 €	1419041	109.70 €	1719041	109.70 €	3219041	170.00 €	1819041	188.00 €	1919041	240.00 €
71,00 mm.●	0818841	75.25 €	1018841	75.25 €	1418841	109.70 €	1718841	109.70 €	3218841	170.00 €	1818841	188.00 €	1918841	240.00 €
63,00 mm.●	0818641	75.25 €	1018641	75.25 €	1418641	109.70 €	1718641	109.70 €	3218641	170.00 €	1818641	188.00 €	1918641	240.00 €
56,00 mm.●	0818441	75.25 €	1018441	75.25 €	1418441	109.70 €	1718441	109.70 €	3218441	170.00 €	1818441	188.00 €	1918441	240.00 €
53,00 mm.●	0818241	75.25 €	1018241	75.25 €	1418241	109.70 €	1718241	109.70 €	3218241	170.00 €	1818241	188.00 €	1918241	240.00 €
50,00 mm.●	0818041	75.25 €	1018041	75.25 €	1418041	109.70 €	1718041	109.70 €	3218041	170.00 €	1818041	188.00 €	1918041	240.00 €
45,00 mm.●	0817841	75.25 €	1017841	75.25 €	1417841	109.70 €	1717841	109.70 €	3217841	156.50 €	1817841	177.00 €	1917841	234.50 €
40,00 mm.●	0817641	75.25 €	1017641	75.25 €	1417641	109.70 €	1717641	109.70 €	3217641	156.50 €	1817641	177.00 €	1917641	234.50 €
37,50 mm.●	0817441	75.25 €	1017441	75.25 €	1417441	109.70 €	1717441	109.70 €	3217441	156.50 €	1817441	177.00 €	1917441	234.50 €
35,50 mm.●	0817241	75.25 €	1017241	75.25 €	1417241	109.70 €	1717241	109.70 €	3217241	156.50 €	1817241	177.00 €	1917241	234.50 €
31,50 mm.●	0817041	75.25 €	1017041	75.25 €	1417041	109.70 €	1717041	109.70 €	3217041	156.50 €	1817041	177.00 €	1917041	234.50 €
28,00 mm.●	0816841	75.25 €	1016841	75.25 €	1416841	109.70 €	1716841	109.70 €	3216841	156.50 €	1816841	177.00 €	1916841	234.50 €
26,50 mm.●	0816641	75.25 €	1016641	75.25 €	1416641	109.70 €	1716641	109.70 €	3216641	156.50 €	1816641	177.00 €	1916641	234.50 €
25,00 mm.●	0816441	75.25 €	1016441	75.25 €	1416441	109.70 €	1716441	109.70 €	3216441	156.50 €	1816441	177.00 €	1916441	234.50 €
22,40 mm.●	0816241	75.25 €	1016241	75.25 €	1416241	109.70 €	1716241	109.70 €	3216241	156.50 €	1816241	177.00 €	1916241	234.50 €
20,00 mm.●	0816041	75.25 €	1016041	75.25 €	1416041	109.70 €	1716041	109.70 €	3216041	156.50 €	1816041	177.00 €	1916041	234.50 €
19,00 mm.●	0815841	75.25 €	1015841	75.25 €	1415841	109.70 €	1715841	109.70 €	3215841	156.50 €	1815841	177.00 €	1915841	234.50 €
18,00 mm.●	0815641	75.25 €	1015641	75.25 €	1415641	109.70 €	1715641	109.70 €	3215641	156.50 €	1815641	177.00 €	1915641	234.50 €
16,00 mm.●	0815441	75.25 €	1015441	75.25 €	1415441	109.70 €	1715441	109.70 €	3215441	156.50 €	1815441	177.00 €	1915441	234.50 €
14,00 mm.●	0815241	75.25 €	1015241	75.25 €	1415241	109.70 €	1715241	109.70 €	3215241	156.50 €	1815241	177.00 €	1915241	234.50 €
13,20 mm.●	0815041	75.25 €	1015041	75.25 €	1415041	109.70 €	1715041	109.70 €	3215041	156.50 €	1815041	177.00 €	1915041	234.50 €
12,50 mm.●	0814841	75.25 €	1014841	75.25 €	1414841	109.70 €	1714841	109.70 €	3214841	156.50 €	1814841	177.00 €	1914841	234.50 €
11,20 mm.●	0814641	75.25 €	1014641	75.25 €	1414641	109.70 €	1714641	109.70 €	3214641	156.50 €	1814641	177.00 €	1914641	234.50 €

**STAINLESS STEEL TEST SIEVES
ACCORDING ISO 3310.2
PERFORATED PLATES - ROUND OPENINGS**

Opening	200 mm.Ø x 50 mm.H		203 mm.Ø x 50 mm.H		300 mm.Ø x 70 mm.H		305 mm.Ø x 80 mm.H		350 mm.Ø x 60 mm.H		400 mm.Ø x 100 mm.H		450 mm.Ø x 100 mm.H	
	M-LB code	Price	M-LB code	Price	M-LB code	Price	M-LB code	Price	M-LB code	Price	M-LB code	Price	M-LB code	Price
10,00 mm.●	0814441	75.25 €	1014441	75.25 €	1414441	109.70 €	1714441	109.70 €	3214441	156.50 €	1814441	177.00 €	1914441	234.50 €
9,50 mm.●	0814241	75.25 €	1014241	75.25 €	1414241	109.70 €	1714241	109.70 €	3214241	156.50 €	1814241	177.00 €	1914241	234.50 €
9,00 mm.●	0814041	75.25 €	1014041	75.25 €	1414041	109.70 €	1714041	109.70 €	3214041	156.50 €	1814041	177.00 €	1914041	234.50 €
8,00 mm.●	0813841	75.25 €	1013841	75.25 €	1413841	109.70 €	1713841	109.70 €	3213841	156.50 €	1813841	177.00 €	1913841	234.50 €
7,10 mm.●	0813641	75.25 €	1013641	75.25 €	1413641	109.70 €	1713641	109.70 €	3213641	156.50 €	1813641	177.00 €	1913641	234.50 €
6,70 mm.●	0813441	75.25 €	1013441	75.25 €	1413441	109.70 €	1713441	109.70 €	3213441	156.50 €	1813441	177.00 €	1913441	234.50 €
6,30 mm.●	0813241	75.25 €	1013241	75.25 €	1413241	109.70 €	1713241	109.70 €	3213241	182.00 €	1813241	203.00 €	1913241	260.50 €
5,60 mm.●	0813041	75.25 €	1013041	75.25 €	1413041	109.70 €	1713041	109.70 €	3213041	182.00 €	1813041	203.00 €	1913041	260.50 €
5,00 mm.●	0812841	75.25 €	1012841	75.25 €	1412841	109.70 €	1712841	109.70 €	3212841	182.00 €	1812841	203.00 €	1912841	260.50 €
4,75 mm.●	0812641	75.25 €	1012641	75.25 €	1412641	109.70 €	1712641	109.70 €	3212641	188.00 €	1812641	208.50 €	1912641	282.00 €
4,50 mm.●	0812441	75.25 €	1012441	75.25 €	1412441	109.70 €	1712441	109.70 €	3212441	188.00 €	1812441	208.50 €	1912441	282.00 €
4,00 mm.●	0812241	75.25 €	1012241	75.25 €	1412241	109.70 €	1712241	109.70 €	3212241	188.00 €	1812241	208.50 €	1912241	282.00 €

SPECIAL ROUND OPENINGS

Opening	M-LB code	Price	M-LB code	Price	M-LB code	Price	M-LB code	Price	M-LB code	Price	M-LB code	Price	M-LB code	Price
32,00 mm.●	0817141	75.25 €	1017141	75.25 €	1417141	109.70 €	1717141	109.70 €	3217141	156.50 €	1817141	177.00 €	1917141	234.50 €
30,00 mm.●	0816941	75.25 €	1016941	75.25 €	1416941	109.70 €	1716941	109.70 €	3216941	156.50 €	1816941	177.00 €	1916941	234.50 €
24,00 mm.●	0816341	75.25 €	1016341	75.25 €	1416341	109.70 €	1716341	109.70 €	3216341	156.50 €	1816341	177.00 €	1916341	234.50 €
22,00 mm.●	0816141	75.25 €	1016141	75.25 €	1416141	109.70 €	1716141	109.70 €	3216141	156.50 €	1816141	177.00 €	1916141	234.50 €
21,00 mm.●	0816081	75.25 €	1016081	75.25 €	1416081	109.70 €	1716081	109.70 €	3216081	156.50 €	1816081	177.00 €	1916081	234.50 €
15,00 mm.●	0815381	75.25 €	1015381	75.25 €	1415381	109.70 €	1715381	109.70 €	3215381	156.50 €	1815381	177.00 €	1915381	234.50 €
14,50 mm.●	0815341	75.25 €	1015341	75.25 €	1415341	109.70 €	1715341	109.70 €	3215341	156.50 €	1815341	177.00 €	1915341	234.50 €
13,00 mm.●	0814941	75.25 €	1014941	75.25 €	1414941	109.70 €	1714941	109.70 €	3214941	156.50 €	1814941	177.00 €	1914941	234.50 €
12,00 mm.●	0814741	75.25 €	1014741	75.25 €	1414741	109.70 €	1714741	109.70 €	3214741	156.50 €	1814741	177.00 €	1914741	234.50 €
11,00 mm.●	0814541	75.25 €	1014541	75.25 €	1414541	109.70 €	1714541	109.70 €	3214541	156.50 €	1814541	177.00 €	1914541	234.50 €

SPECIAL ROUND OPENINGS

Opening	200 mm.Ø x 50 mm.H		203 mm.Ø x 50 mm.H		300 mm.Ø x 70 mm.H		305 mm.Ø x 80 mm.H		350 mm.Ø x 60 mm.H		400 mm.Ø x 100 mm.H		450 mm.Ø x 100 mm.H	
	M-LB code	Price	M-LB code	Price	M-LB code	Price	M-LB code	Price	M-LB code	Price	M-LB code	Price	M-LB code	Price
10,20 mm.●	0814471	75.25 €	1014471	75.25 €	1414471	109.70 €	1714471	109.70 €	3214471	156.50 €	1814471	177.00 €	1914471	234.50 €
9,30 mm.●	0814141	75.25 €	1014141	75.25 €	1414141	109.70 €	1714141	109.70 €	3214141	156.50 €	1814141	177.00 €	1914141	234.50 €
8,50 mm.●	0813941	75.25 €	1013941	75.25 €	1413941	109.70 €	1713941	109.70 €	3213941	156.50 €	1813941	177.00 €	1913941	234.50 €
7,50 mm.●	0813741	75.25 €	1013741	75.25 €	1413741	109.70 €	1713741	109.70 €	3213741	156.50 €	1813741	177.00 €	1913741	234.50 €
7,00 mm.●	0813541	75.25 €	1013541	75.25 €	1413541	109.70 €	1713541	109.70 €	3213541	156.50 €	1813541	177.00 €	1913541	234.50 €
6,50 mm.●	0813341	75.25 €	1013341	75.25 €	1413341	109.70 €	1713341	109.70 €	3213341	156.50 €	1813341	177.00 €	1913341	234.50 €
6,00 mm.●	0813141	75.25 €	1013141	75.25 €	1413141	109.70 €	1713141	109.70 €	3213141	156.50 €	1813141	177.00 €	1913141	234.50 €
3,55 mm.●	0812041	97.70 €	1012041	97.70 €	1412041	165.00 €	1712041	165.00 €	3212041	230.00 €	1812041	265.00 €	1912041	307.00 €
3,40 mm.●	0811941	97.70 €	1011941	97.70 €	1411941	165.00 €	1711941	165.00 €	3211941	230.00 €	1811941	265.00 €	1911941	307.00 €
3,15 mm.●	0811641	97.70 €	1011641	97.70 €	1411641	165.00 €	1711641	165.00 €	3211641	230.00 €	1811641	265.00 €	1911641	307.00 €
3,00 mm.●	0811541	97.70 €	1011541	97.70 €	1411541	165.00 €	1711541	165.00 €	3211541	230.00 €	1811541	265.00 €	1911541	307.00 €
2,95 mm.●	0811471	97.70 €	1011471	97.70 €	1411471	165.00 €	1711471	165.00 €	3211471	230.00 €	1811471	265.00 €	1911471	307.00 €
2,80 mm.●	0811441	97.70 €	1011441	97.70 €	1411441	165.00 €	1711441	165.00 €	3211441	230.00 €	1811441	265.00 €	1911441	307.00 €
2,65 mm.●	0811341	97.70 €	1011341	97.70 €	1411341	165.00 €	1711341	165.00 €	3211341	230.00 €	1811341	265.00 €	1911341	307.00 €
2,50 mm.●	0811241	97.70 €	1011241	97.70 €	1411241	165.00 €	1711241	165.00 €	3211241	230.00 €	1811241	265.00 €	1911241	307.00 €
2,35 mm.●	0811041	97.70 €	1011041	97.70 €	1411041	165.00 €	1711041	165.00 €	3211041	230.00 €	1811041	265.00 €	1911041	307.00 €
2,24 mm.●	0810841	97.70 €	1011041	97.70 €	1411441	165.00 €	1711741	165.00 €	3213241	230.00 €	1811841	265.00 €	1911941	307.00 €
2,00 mm.●	0810641	97.70 €	1010641	97.70 €	1410641	165.00 €	1710641	165.00 €	3210641	230.00 €	1810641	265.00 €	1910641	307.00 €
1,80 mm.●	0810441	158.00 €	1010441	158.00 €	1410441	195.00 €	1710441	195.00 €	3210441	230.00 €	1810441	265.00 €	1910441	307.00 €
1,60 mm.●	0810041	158.00 €	1010041	158.00 €	1410041	195.00 €	1710041	195.00 €	3210041	230.00 €	1810041	265.00 €	1910041	307.00 €
1,50 mm.●	0809941	158.00 €	1009941	158.00 €	1409941	195.00 €	1709941	195.00 €	3209941	230.00 €	1809941	265.00 €	1909941	307.00 €
1,40 mm.●	0809841	158.00 €	1009841	158.00 €	1409841	195.00 €	1709841	195.00 €	3209841	230.00 €	1809841	265.00 €	1909841	307.00 €
1,25 mm.●	0809641	158.00 €	1009641	158.00 €	1409641	195.00 €	1709641	195.00 €	3209641	230.00 €	1809641	265.00 €	1909641	307.00 €
1,00 mm.●	0809041	158.00 €	1009041	158.00 €	1409041	195.00 €	1709041	195.00 €	3209041	230.00 €	1809041	265.00 €	1909041	307.00 €
500 µm.●	0807441	158.00 €	1007441	158.00 €	1407441	195.00 €	1707441	195.00 €	3207441	230.00 €	1807441	265.00 €	1907441	307.00 €

STAINLESS STEEL TEST SIEVES
DIAMETER : 200 mm.Ø
HIGHT : 50 mm.H
STAINLESS STEEL PERFORATED PLATE
ELONGATED PERFORATION

Opening	M-LB code	Price	Standard
1,00 x 20,00 mm	0830081	91.50 €	ISO 5223
1,20 x 20,00 mm	0832781	91.50 €	
1,40 x 20,00 mm	0833081	91.50 €	
1,50 x 20,00 mm	0830181	91.50 €	
1,60 x 20,00 mm	0832881	91.50 €	
1,70 x 20,00 mm	0830281	91.50 €	ISO 5223
1,75 x 20,00 mm	0833181	91.50 €	
1,80 x 20,00 mm	0830381	91.50 €	ISO 5223
1,90 x 20,00 mm	0830481	91.50 €	ISO 5223
2,00 x 20,00 mm	0830581	91.50 €	ISO 5223
2,10 x 20,00 mm	0830681	91.50 €	
2,20 x 20,00 mm	0830781	91.50 €	ISO 5223
2,25 x 20,00 mm	0830881	91.50 €	
2,30 x 20,00 mm	0830981	91.50 €	
2,40 x 20,00 mm	0831081	91.50 €	
2,50 x 20,00 mm	0831181	91.50 €	ISO 5223
2,60 x 20,00 mm	0831281	91.50 €	
2,75 x 20,00 mm	0831381	91.50 €	
2,80 x 20,00 mm	0831481	91.50 €	ISO 5223
2,90 x 20,00 mm	0831581	91.50 €	
3,00 x 20,00 mm	0831681	91.50 €	
3,20 x 20,00 mm	0831781	91.50 €	
3,25 x 20,00 mm	0833281	91.50 €	
3,55 x 20,00 mm	0831881	91.50 €	ISO 5223
3,75 x 20,00 mm	0832981	91.50 €	
4,00 x 20,00 mm	0832581	91.50 €	
4,50 x 25,00 mm	0832181	91.50 €	
4,76 x 25,00 mm	0831981	91.50 €	
5,00 x 30,00 mm	0832281	91.50 €	
5,50 x 30,00 mm	0832681	91.50 €	
5,60 x 25,00 mm	0832081	91.50 €	

STAINLESS STEEL TEST SIEVES
DIAMETER : 200 mm.Ø
HIGHT : 50 mm.H
NYLON MESH

Opening	M-LB Code	Price
5µm.#	0800211	124.00 €
10µm.#	0800411	102.40 €
20µm.#	0800611	102.40 €
20µm.#	0800631	102.40 €
25µm.#	0800811	102.40 €
32µm.#	0801011	102.40 €
36µm.#	0801211	102.40 €
38µm.#	0801411	102.40 €
40µm.#	0801611	102.40 €
45µm.#	0801811	102.40 €
50µm.#	0802011	102.40 €
50µm.#	0802031	102.40 €
53µm.#	0802211	102.40 €
56µm.#	0802411	102.40 €
63µm.#	0802611	102.40 €
71µm.#	0802811	102.40 €
75µm.#	0803011	102.40 €
80µm.#	0803211	102.40 €
90µm.#	0803411	102.40 €
100µm.#	0803611	102.40 €
106µm.#	0803811	79.00 €
112µm.#	0804011	79.00 €
125µm.#	0804211	79.00 €
140µm.#	0804411	79.00 €
150µm.#	0804611	79.00 €
160µm.#	0804811	79.00 €
180µm.#	0805011	79.00 €
200µm.#	0805211	79.00 €
200µm.#	0805231	79.00 €
212µm.#	0805411	79.00 €
224µm.#	0805611	79.00 €
250µm.#	0805811	79.00 €
280µm.#	0806011	79.00 €
300µm.#	0806211	79.00 €
315µm.#	0806411	79.00 €
355µm.#	0806611	79.00 €
400µm.#	0806811	79.00 €
425µm.#	0807011	79.00 €
450µm.#	0807211	79.00 €
500µm.#	0807411	79.00 €
560µm.#	0807611	79.00 €
600µm.#	0807811	79.00 €
630µm.#	0808011	79.00 €
710µm.#	0808211	79.00 €
800µm.#	0808411	79.00 €
850µm.#	0808611	79.00 €
900µm.#	0808811	79.00 €
1,00mm.#	0809011	79.00 €
2,00mm.#	0810631	79.00 €